

Otto Bock®

QUALITY FOR LIFE



C-Leg® Microprocessor Knee

Technology that Transforms Lives

Since it was introduced in the U.S. in 1999, more than 30,000 people worldwide have taken advantage of the unparalleled stability offered by the C-Leg® microprocessor-controlled knee—making this technology the Standard of Care.

With the C-Leg® you can speed up, slow down, take on hills and go down stairs step-over-step—all with a secure, more natural gait.



Features and Benefits

- **Transformative Technology.** The revolutionary microprocessor that controls the C-Leg® receives feedback from multiple sensors, allowing the knee to anticipate your next move and make adjustments in real-time.
- **Always Thinking About Your Security—So You Don't Have To.** The C-Leg® offers unprecedented stability—and flexibility. You can speed up, slow down, and go down stairs step-over-step with the assurance that the C-Leg® is there for you.
- **Stumble Recovery.** Whenever the C-Leg® senses that you're in an insecure position—such as when you trip or stumble—it will stiffen to provide the support you need to recover.
- **Multiple Modes for Multiple Uses.** The C-Leg® offers several functional modes. The first two modes are set by your prosthetist for walking and other favorite activities, such as driving or riding a bike.
- **Stand at Ease.** A new Static Stance mode lets you lock the leg yourself between 7° and 70°, helping you use less energy no matter what you stand for! A vibration alert lets you know when Static Stance is engaged and disengaged.
- **Switch Over.** Switching between modes is as easy as bouncing on your toe—or a press of the wireless remote control.

Is the C-Leg® Technology Right for Me?

The C-Leg® is designed for more active users who can take advantage of its Dynamic Swing Control and Second Mode, while our companion knee joint, the C-Leg® Compact, is right for those who simply need the best in stability. Both knees can accommodate patients who weigh up to 275 lbs/125 kg.

If you have any questions about whether C-Leg® technology is right for you, please talk to your prosthetist. Don't have a prosthetist? Please visit our website to find a prosthetist in your area, www.ottobockus.com/cleg.

Ordering Information

For more information call Otto Bock HealthCare at 800.328.4058 (U.S.) or 800.665.3327 (Canada). You can also find more information online at www.ottobockus.com (U.S.) or www.ottobock.ca (Canada).

Warranty and Service

36-month warranty. Service 24 months after purchase. Please refer to the Instructions for Use for complete information.

